

Battery charger upgrade for Fortress 1700

Starting September 2011, we will begin delivering new Fortress 1700 scooters with an upgraded 8 AMP charger in replacement of our current 5 AMP charger.

Upcoming Events:

Rehab Equipment Expo
Seaforth Armouries
Vancouver, BC
Sep 13, 2011

**Health Care Innovations
Toronto**
Toronto Congress Centre
Toronto, ON
Nov 14 - 15, 2011
<http://toronto.healthcareinnovations.ca>

**Health Care Innovations
Calgary**
Calgary Telus Convention Centre
Calgary, AB
Nov 29, 2011
<http://calgary.healthcareinnovations.ca>

**28th International Seating
Symposium**
Westin Bayshore, Vancouver, BC
Mar 6-9, 2012
<http://www.interprofessional.ubc.ca/iss>

**Handicare Educational
Workshop - Barrie, ON**
Liberty North
Oct 5, 2011
RSVP REQUIRED, please see page 2 of
newsletter for details

**Handicare Educational
Workshop - London, ON**
Best Western Stoneridge Inn &
Conference Centre
Oct 13, 2011
RSVP REQUIRED, please see page 2 of
newsletter for details

FORTRESS WINNER SCOOTER UPDATE:

We would like to take a moment to inform you about an update to the Fortress Winner scooter. All Winner scooters will now come with a new transaxle. This new transaxle performs extremely well with our 140 amp from PG Drives Technology. The motor itself an extra-large 4-pole design that continues to offer a 15km/h top speed just like the previous transaxle offered.



Additional details:

- All serial numbers that begin with "WP" and onward will contain the new transaxle
- The "Alpine Drive System" for the Fortress Winner remains un-changed - only the standard Winner transaxle is new
- We will still be able to provide service parts for the older generation transaxle
- The freewheel lever on the new transaxle operates in the opposite direction as the previous generation transaxle - as such we have changed the freewheel label accordingly

Parts Information:

As a result of this transaxle change there are also a few part number changes. These are outlined in the table below:

Description	New Part #
Transaxle including mounting accessories Winner S-Drive	
Transaxle Winner std. S-Drive CT	1012904
EM Brake set CT	1011584
Motor brushes 17.9x9.9 CT	1011585
Freewheel handle Winner CT	1011586
Key Winner CT	1012123
M10 x 100 8.8 VZ DIN 931	1012699
Transaxle heavy duty including mounting accessories Winner S-Drive	Unchange

Upcoming Educational Workshop:

Building Your Case: When and how to justify the use of Type 3 Manual Wheelchairs or Type 5 Tilt /Dynamic Tilt Chairs -All aspects from clinical application to funding. By Sheila Buck, B.Sc.(OT), Reg. Ont., ATP

Sheila Buck is a well-known Occupational Therapist, educator, and now author, who has been working in the Seating & Mobility field for over 20 years. Sheila will take the mystery out of prescribing Category 3 Rehab Chairs and Category 5 Tilt Chairs by reviewing the benefits of each type of chair related to client presentation and diagnosis.

Barrie, ON

Oct 5, 2011

Time: 12:00pm – 4:00pm

London, ON

Oct 13, 2011

Time: 12:00pm – 4:00pm

If you're interested in any of the above workshops, please RSVP through
Tel: 1-800-514-7061 or
Email: info@handicare.ca

A light lunch will be provided

WE ARE NOW ONLINE!

follow us on
twitter



@HandicareCanada

and become our



FACEBOOK fan on the

Handicare Canada fan page!

SHARE YOUR STORY

Handicare loves hearing from the community. Whether it is an assessment experience, success story, end user case study, challenges, we would love to hear from you.

Please email your stories to: crystal.wong@handicare.ca



Looking for someone to CALL or to schedule an IN-SERVICE? Please see below for sales staff contacts:

STAFF	DEPARTMENT	EXT#	TEL#
Frank Malta	Sales - National Sales Manager	228	416-998-8012
John Russell	Sales - British Columbia	304	778-837-1197
Nigel Moore	Sales - South Western Ontario +GTA	303	416-706-0204
Wayne Parent	Sales - Northern Ontario	305	705-627-5329
Robert Price	Sales - Eastern Ontario + GTA	306	416-709-3880
Steve Chipun	Sales - Rehab Mobility Specialist	236	416-271-9677
Frank Agozzino	Sales - Quebec	302	514-238-6696
Marc Gagnon	Sales - Quebec - Rehab	301	418-930-2053
Indra Suratt	Customer Service	242	416-739-8333
Hala Andary	Customer Service (English/French)	221	416-739-8333