






# Handicare 1000 Perch

Straight Stairlifts

# Handicare 1000 Perch

## Straight Stairlifts



| Perch seat                             |                                      | Dimensions    |
|--|--------------------------------------|---------------|
| A                                      | Minimum swivel radius*               | -             |
| B                                      | Wall to inner seat back              | 2.95in        |
| C                                      | Seat depth                           | 12.59in       |
| D                                      | Minimum folded width to footplate    | 15.15in       |
| E                                      | Minimum open width to edge footplate | 25.59in       |
| F                                      | Armrest width – external             | 21.65in       |
| G                                      | Armrest width – internal             | 18.30in       |
| H                                      | Seat back height                     | 17.91in       |
| I                                      | Footplate to seat height             | 26.57-34.44in |
| J                                      | Minimum footplate height             | 2.75in        |
| Minimum track intrusion into staircase |                                      | 4.52in        |
| Maximum weight limit                   |                                      | 300lbs        |
| Speed                                  |                                      | 0.07m/sec     |
| Staircase angle                        |                                      | 27 - 45°      |

\*Fixed seat

Some of the above dimensions depend on the characteristics of the staircase in the home. Please consult the sales manual for further information about these dimensions and features of this seat. All information is provided in good faith but does not constitute a guarantee. Handicare reserves the right to alter designs and specifications without prior notice.

**HANDICARE**

10888 Metro Court St. Louis, MO 63043 | 866.891.6502 | [www.handicareusa.com](http://www.handicareusa.com)  
81 Romina Drive Concord, ON L4K 4Z9 | 888.637.8155 | [www.handicare.ca](http://www.handicare.ca)