















**LimbSling**  
SystemRoMedic  
Disposable  
REF: 49292000  
Material: Non woven/  
polypropylene/  
Band: Polyester  
SWL: 150kg/330lbs  
CE

Visual inspection  
Visual inspection  
Read the manual  
Read the manual

**DO NOT WASH!**  
handicare  
Handicare AB  
Box 10000  
S-201 08 Stockholm  
Sweden  
www.handicare.com

## Disposable - but not for single use

All of our disposable transfer aids may be used for one user throughout their hospital stay. Handicare's disposable transfer aids are safe, cost effective and convenient options facilitating infection control efforts.

Our disposable transfer aids include the following items and others on their way!



## Wipeable transfer aids in several materials

Our wipeable transfer aids are available in several materials and configurations.



For updates and more information, please visit our website, [www.handicare.com](http://www.handicare.com).

Handicare provides a complete solution for easy transfers. The concept encompasses four product areas that meet most needs in different transfer situations.

**Transfer**, between two locations. **Positioning**, within the same location.

**Support**, for mobility. **Lifting**, both manual and mechanical.

.FUSPPVSU  
4U-PVT.0

.MMXBW6  
PDPSE0/-,7

XXXIBOEDBSFVTBDPN

XXXIBEDBSDB