

SafeHandlingSheet Dimensions

Name	SafeHandlingSheet				
Part Number	47522007	47542007	47542007	47502009	547810
Description	Polyester net, L	Polyester net/ WendyLett4Way, L	Polyester, XL	Polyester, reinforced, XXL	Disposable, one size
Unit	1 piece	1 piece	1 piece	1 piece	10 pieces
Width	45.3in	45.3in	45.3in	51.2in	52in
Length	94.5in	94.5in	94.5in	94.5in	78in
Recommended User Weight	110lbs – 660lbs	110lbs – 660lbs	110lbs – 880lbs	110lbs – 1257lbs	110lbs – 1000lbs
Safe Working Load	660lbs	660lbs	880lbs	1257lbs	1000lbs