

WendyLett Dimensions

Name	WendyLett						
Part Number	1609	1614	1629	1634	1639	1667	1668
Description	White		Gray/white striped				
Unit	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece
Length	39.4"	55.1"	39.4"	55.1"	78.7"	39.4"	55.1"
Width	78.7"	78.7"	78.7"	78.7"	78.7"	86.6"	86.6"
Sliding Surface	23.6"	23.6"	23.6"	23.6"	27.6"	39.4"	39.4"

Name	WendyLett					
Part Number	1669	1641	1651	1661	1671	1681*
Description	Gray/white striped	Fitted Sheet w/ Elastic, gray/white striped				
Unit	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece
Length	78.7"	78.7"	78.7"	78.7"	78.7"	78.7"
Width	86.6"	35.4"	41.3"	47.2"	55.1"	145.7"
Sliding Surface	39.4"	27.6"	31.5"	37.4"	43.3"	27.6"

*extra wide to allow for optional placing of sliding surface

Name	WendyLett2Way			Accessories
Part Number	1620	1631	1636	70200007
Description	Draw sheet, white		Draw sheet, gray/white	RoMedicClips
Unit	1 piece	1 piece	1 piece	2 pieces
Length	35.4"	39.4"	55.1"	-
Width	70.9"	78.7"	78.7"	-

Name	WendyLett4Way					
Part Number	1637	1647	1648	1649	1657*	1652
Description	Sliding draw sheet, gray/white checked		Sliding draw sheet/ incontinence protection with extra fabric on each side, gray/white checked	Sliding draw sheet/ incontinence protection with 7 handles on each side, gray	Sliding draw sheet with 7 handles on each side, gray/white checked	Bariatric use, gray/white checked
Unit	1 piece	1 piece	1 piece	1 piece	1 piece	1 piece
Length	55.1"	78.7"	29.5"	55.1"	55.1"	55.1"
Width	78.7"	78.7"	78.7"	47.2"	74.8"	102.4"

*Can be used in combination with mobile or stationary lifts