



# EVE! Bath and Shower Lift

Modern Design Meets Safety and Versatility



## Safe and Comfortable Shower and Bathing Experiences

The EVE! Bath and Shower Lift features an ergonomic design allowing caregivers to provide safe and comfortable shower and bathing experiences for residents. The unique design of the EVE! makes its extremely versatile and easy to use with simple integration in any existing care bathtub. The padded seat with hygiene opening provides the resident with ultimate comfort while the dual controls on the pillar and handset allow for easy one or two caregiver operation. Additional safety features include an integrated service module reporting any technical problems and required maintenance.

Battery and charger are backwards compatible across all Handicare by Beka mobile lifts.

# Caregiver and Resident Benefits



Safety steering handles enclose the resident ensuring safe and easy maneuvering.



Cable remote control with display of the charging condition of the battery, failure diagnostics, service.



Comfortable shaped seat with hygiene opening.



Built-in door with emergency stop allows for easy battery access and removal.



Closed, easy-to-maneuver double castors. Both back wheels include braking system.



Lowers to a height allowing the resident to sit down on the seat with ease.



Integrated digital scale powered by the main battery (optional).



The EVE! fits in any care bathtub.



Aluminum components make the lift extremely lightweight.

# EVE! Bath and Shower Lift Specifications

## Part Numbers

### Safe Working Load (SWL)

185 kg (408 lb)

### Dimensions (L x W x H)

Length: 950 mm (37.4 in)

Width: 680 mm (26.8 in)

Height: 1240 mm (48.8 in)

Height, seat: 495 mm (19.5 in), minimum

Height, chassis: 140 mm (5.5 in)

Free space under carrier frame: 30 mm (1.18 in)

### Lifting Range

495 – 1065 mm (19.5 – 41.9 in)

### Castors

Two double castors, front: 100 mm dia. (3.9 in)

Two double castors, rear, fixed: 100 mm dia. (3.9 in)

### Power Supply

Removable 24V battery

### Protection Class

IPX4

### Product Weight Without Packaging

63 kg (139 lb)

For additional technical data, please refer to the user manual.

**920107000:** EVE! Bathing Lift

**920108000:** EVE! Bathing Lift with Digital Scale

## Accessories / Options

Description	Part No.
Digital scale (retrofit)	920107500
Replacement battery	920107602
External wall charging station	920107702



### HANDICARE

485 Millway Avenue, Unit 2, Concord, ON, L4K 3V4

Toll Free: 877.304.5438 | Fax: 416.260.0256

[customerservice.canada@handicare.com](mailto:customerservice.canada@handicare.com)

[www.handicare.ca](http://www.handicare.ca)

10888 Metro Court, St. Louis, MO 63043

Toll Free: 866.891.6502 | Fax: 314.692.7858

[customerservice.northamerica@handicare.com](mailto:customerservice.northamerica@handicare.com)

[www.handicareusa.com](http://www.handicareusa.com)

Are you ready to invest in innovation  
and simplify bathing for your residents?

To learn more, call 1-888-637-8155 or  
visit [www.handicare.ca](http://www.handicare.ca) / [handicareusa.com](http://handicareusa.com)